

Reconditioned North American Drager 2C



Shown with (3) 19.1 style vaporizers.

Manufactured from 1991-2001 the Narkomed 2C features.

O₂ and N₂O standard

Monitoring: Oxygen(inspired), Volume & Pressure

ORC* (Hypoxic Guard)

Breathing circuit: Traditional size and configuration of absorber head with higher-volume breathing circuit. Accepts canister or 10 liter absorber. May be used with non-rebreathing circuits.

Scavenger: Open or Closed

Ventilator: AVE2+ ventilator w/ ascending bellows standard-Modes volume control (VCV), Manual/Spontaneous. Mechanical "bag/vent" switch.

E-Cylinder and central gas connections standard

Flowmeters: Traditional needle valves with glass flow tube display

Vaporizers:Removable variable By Pass Vaporizers

Physical Specs:

Drawers: 2 and/or 3 drawer

- 5" shorter than 2B

Displays: Removable flat plate displays mounted to the articulating arm: Displays volume, pressure, and oxygen data through trend alarm system.