

## Reconditioned North American Drager Narkomed GS



The Narkomed GS was manufactured from 1997 to 2006. The GS is very highly demanded due to its familiar Narkomed design, dependability, and economic value. A dependable and familiar flow head with visible bellows provide reliable performance with Volume, Pressure and Oxygen monitoring.

### Features Include:

- Monitoring Volume, Pressure, inspired Oxygen (VPO), Galvanic cell oxygen analysis, Electronic Wright respirometer style flow sensor
- Ventilator: Ascending pneumatic bellows AV2+ style ventilator without tidal volume compensation. Modes volume control (VCV), Manual/Spontaneous. Mechanical "bag/vent" switch.
- Flowmeters: Traditional needle valve with glass flowtube display. Vaporizers: Variable-bypass vaporizers are used. Vapor 19 or Tec 6. Available in 2 or 3 vaporizer configurations.
- Breathing circuit: Traditional absorber head. The breathing circuit is higher volume (dual absorbent canisters). May use loose granules or prefilled canisters. May be used with non-rebreathing circuits.
- Scavenger: Open, closed, or passive scavenger interface.
- Electrical power failure: Most functions except patient monitoring are preserved during power failure, as this is a basic pneumatic-mechanical machine. The VPO display will be lost after battery exhaustion. GS offers a large 3 drawer cabinet, writing tray and auxiliary flow meter. Owners of the Narkomed GS must supply patient monitors & gas analysis.